

Blauman[®]Effect

The **laser revolution** that removes imperfections
and enhances the skin's radiance



M milesman
blauman



100% *Manufactured
in Europe*

Attain the **BlaumanEffect**[®] on your skin

These are the results claimed by professionals who are already using the **Blueman blue laser**.

After revitalising the skin with the laser they observed:

- ✓ Increased radiance on the skin.
- ✓ More hydrated face.
- ✓ Smoother face.
- ✓ Tighter tone.
- ✓ Skin more receptive to receiving other treatments.
- ✓ Cleaner pores.
- ✓ Removal of microlesions such as milia, sebaceous hyperplasias, telangiectasias.



The ultimate **acne treatment**



In a single session, the Milesman Blauman blue laser reduces 69% of active acne lesions in just 72 hours and 83% in 7 days.

Study conducted by
Milesman I Gaias Research 2020



Scan

Removes imperfections in seconds

Blauman blue laser is the perfect device to remove all your customers' small, benign, aesthetic lesions. Warts, ruby spots or spider veins disappear in very few sessions and there is no need for a scalpel or anesthesia.

How it works

Thanks to the laser equipment we can vaporize the tissue treated (excess skin, skin lesions, stretch marks, scars), and instantly remove it without transferring heat to the adjacent tissues.

The kind of lesions we can treat

Vascular lesions

There are two main kinds of vascular lesions that we can treat with our Blauman blue laser:



SPIDER VEINS

Pigmented lesions

They are very common and quickly and safely removed with the Blauman blue laser:



NEVI

Other lesions

With the Blauman blue laser we can also remove other kinds of skin lesions:



ACROCHORDONS



RUBY SPOTS



SEBORRHEIC KERATOSIS



LENTIGINES

Technical specifications

Wavelength

450 nm

Operating mode

Pulsed

Spot size

Two interchangeable heads
0.3 mm and 1 cm²

Maximum power

10 W

Pulse width

20 to 200 ms

Laser dimensions

25 cm x 5 cm

Laser weight

2.2 kg



Benefits of the blue laser



Skin assessments



Kills bacteria

Has an antiseptic effect on the treated area.



No anesthetic



No bleeding

The laser photocoagulates the treated area, completely eliminating the chance of bleeding.



Fast

It takes less than 10 minutes, on average, to remove one lesion in one session.



Precise

The 0.3 mm tip helps treat the precise area affected.



Portable

2.2 kg in weight and 25 cm in diameter.